



**Product description**

Order number	83.3912.300
Product description	Cell culture flask, T-175, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, 2-position screw cap, TC Tested, 5 piece(s)/bag

**Product characteristics**

Cap	2-position screw cap
Cell type	challenging adherent cells
Growth area	175 cm <sup>2</sup>
Volume of work	50 ml - 125 ml
Vessel type	flask
Coding colour	yellow

**Material & colours**

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	40
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	40
Piece(s) / pallet	960
Depth of case	522 mm
Width of case	267 mm
Height of case	450 mm
Case volume	0.0627 cbm
Weight of product	0.1004 kg
Weight of case	5.56 kg
EAN of inner box	4038917337533
EAN case	4038917337526

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

