



Product description

Order number	83.3920.005
Product description	Cell culture plate, 6 well, material: PS, surface: Standard, for Adherent cells, coding colour: red, flat base, TC Tested, 5 piece(s)/bag

Product characteristics

Label/ Print	with plastic label
Closure type	slip-on lid
Cell type	Adherent cells
Base shape	flat base
Number of wells	6
Growth area	9.07 cm ²
Volume of work	4 ml
Vessel type	Plate
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Cell culture plate, 6 well, surface: Standard, flat base

Page 2

Size

Max. volume	15.75 ml
Width of product	85.2 mm
Height of product	22.4 mm
Length of product	127.8 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polystyrene (PS)

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	2000
Depth of case	567 mm
Width of case	372 mm
Height of case	342 mm
Case volume	0.0721 cbm
Weight of product	0.0449 kg
Weight of case	7.26 kg
EAN of inner box	4038917360135
EAN case	4038917360128

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-03-11

This document was prepared by EDP support and is valid without signature.



SARSTEDT