



Product description

| | |
|---------------------|---|
| Order number | 83.3920.500 |
| Product description | Cell culture plate, 6 well, material: PS, surface: Suspension, for suspension cells, coding colour: green, flat base, TC Tested, 1 piece(s)/blister |

Product characteristics

| | |
|-----------------|----------------------|
| Label/ Print | with plastic label |
| Closure type | slip-on lid |
| Cell type | suspension cells |
| Base shape | flat base |
| Number of wells | 6 |
| Growth area | 9.07 cm ² |
| Volume of work | 4 ml |
| Vessel type | Plate |
| Coding colour | green |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|----------|
| Max. volume | 15.75 ml |
| Width of product | 85.2 mm |
| Height of product | 22.4 mm |
| Length of product | 127.8 mm |

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |
| Closure material | Polystyrene (PS) |

Purity & certification

| | |
|------------------|-------------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Purity standard | TC Tested |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 50 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 50 |
| Piece(s) / pallet | 2400 |
| Depth of case | 332 mm |
| Width of case | 297 mm |
| Height of case | 275 mm |
| Case volume | 0.0271 cbm |
| Weight of product | 0.0451 kg |
| Weight of case | 3 kg |
| EAN of inner box | 4038917338400 |
| EAN case | 4038917338394 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

