



## Product description

Order number	83.3920
Product description	Cell culture plate, 6 well, material: PS, surface: Standard, for Adherent cells, coding colour: red, flat base, TC Tested, 1 piece(s)/blister

## Product characteristics

Label/ Print	with plastic label
Closure type	slip-on lid
Cell type	Adherent cells
Base shape	flat base
Number of wells	6
Growth area	9.07 cm <sup>2</sup>
Volume of work	4 ml
Vessel type	Plate
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Cell culture plate, 6 well, surface: Standard, flat base**

**Page 2**

**Size**

Max. volume	15.75 ml
Width of product	85.2 mm
Height of product	22.4 mm
Length of product	127.8 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polystyrene (PS)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	2400
Depth of case	332 mm
Width of case	297 mm
Height of case	275 mm
Case volume	0.0271 cbm
Weight of product	0.0451 kg
Weight of case	3 kg
EAN of inner box	4038917338349
EAN case	4038917338332

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

