



### **Product description**

Order number	83.3921.500
Product description	Cell culture plate, 12 well, material: PS, surface: Suspension, for suspension cells, coding colour: green, flat base, TC Tested, 1 piece(s)/blister

### **Product characteristics**

Label/ Print	with plastic label
Closure type	slip-on lid
Cell type	suspension cells
Base shape	flat base
Number of wells	12
Growth area	3.65 cm <sup>2</sup>
Volume of work	2 ml
Vessel type	Plate
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



### **Size**

Max. volume	6.64 ml
Width of product	85.2 mm
Height of product	22.5 mm
Length of product	127.8 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polystyrene (PS)

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	2400
Depth of case	332 mm
Width of case	297 mm
Height of case	275 mm
Case volume	0.0271 cbm
Weight of product	0.0501 kg
Weight of case	3.25 kg
EAN of inner box	4038917338493
EAN case	4038917338486

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

