



Product description

Order number	83.3922.500
Product description	Cell culture plate, 24 well, material: PS, surface: Suspension, for suspension cells, coding colour: green, flat base, TC Tested, 1 piece(s)/blister

Product characteristics

Label/ Print	with plastic label
Closure type	slip-on lid
Cell type	suspension cells
Base shape	flat base
Number of wells	24
Growth area	1.82 cm ²
Volume of work	1 ml
Vessel type	Plate
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Cell culture plate, 24 well, surface: Suspension,
flat base**

Page 2

Size

Max. volume	3.39 ml
Width of product	85.2 mm
Height of product	22.45 mm
Length of product	127.8 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	2400
Depth of case	332 mm
Width of case	297 mm
Height of case	275 mm
Case volume	0.0271 cbm
Weight of product	0.0569 kg
Weight of case	3.59 kg
EAN of inner box	4038917338585
EAN case	4038917338578

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT