



Product description

Order number	83.3931.101
Product description	TC insert, for 12-well plates, membrane: PET, transparent, pore size: 1 µm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Suitable for	12-well plates
Base shape	flat base
Growth area	1.1 cm²
Volume of work	0,2 ml - 0,8 ml
Vessel type	TC insert
Pore size	1

Size

Height of product	16.38 mm
-------------------	----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Membrane material	Polyethylene therephthalate (PET)
Frame material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	X-Ray
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	48
Type of smallest subpackaging	blister
Piece(s) / inner box	48
Piece (s) / outer case	1440
Piece(s) / pallet	23040
Depth of box	196 mm
Width of box	104 mm
Height of box	131 mm
Depth of case	410 mm
Width of case	435 mm
Height of case	550 mm
Case volume	0.0981 cbm
Weight of product	0.0037 kg
Weight of case	8.56 kg
EAN of inner box	4038917371193
EAN case	4038917371186

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

