



### Product description

Order number	83.3932.041
Product description	TC insert, for 24-well plates, membrane: PET, transparent, pore size: 0.4 µm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

### Product characteristics

Suitable for	24-well plates
Base shape	flat base
Growth area	0.3 cm <sup>2</sup>
Volume of work	0,1 ml - 0,4 ml
Vessel type	TC insert
Pore size	0.4 µm

### Size

Height of product	16.32 mm
-------------------	----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Membrane material	Polyethylene therephthalate (PET)
Frame material	Polystyrene (PS)

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	X-Ray
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	48
Type of smallest subpackaging	blister
Piece(s) / inner box	48
Piece (s) / outer case	1440
Piece(s) / pallet	23040
Depth of box	196 mm
Width of box	104 mm
Height of box	131 mm
Depth of case	550 mm
Width of case	415 mm
Height of case	440 mm
Case volume	0.1004 cbm
Weight of product	0.0032 kg
Weight of case	7.87 kg
EAN of inner box	4038917371346
EAN case	4038917371339

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

