



Product description

| | |
|---------------------|---|
| Order number | 83.9923.929 |
| Product description | Cell culture tubes, 10 ml, (LxØ): 97 x 16 mm, skirted conical base, PS, TC-treated, for Suspension and mono-layer cultures, sterile, 5 piece(s)/bag |

Product characteristics

| | |
|---------------------------|------------------------------------|
| Graduation | no |
| Cap | standard screw cap |
| Cell type | Suspension and mono-layer cultures |
| Base shape | skirted conical base |
| Volume of work | 10 ml |
| Vessel type | Tube |
| Type of surface treatment | TC-treated |
| Coding colour | red |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Cell culture tubes, (LxØ): 97 x 16 mm, skirted
conical base, TC-treated**

Page 2

Size

| | |
|-------------------|-------|
| Diameter | 16 mm |
| Length of product | 97 mm |

Material & colours

| | |
|-------------------|-------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |
| Closure material | Polyethylene (PE) |
| Colour of cap | red |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 5 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 20000 |
| Depth of case | 428 mm |
| Width of case | 438 mm |
| Height of case | 372 mm |
| Case volume | 0.0697 cbm |
| Weight of product | 0.0084 kg |
| Weight of case | 9.2 kg |
| EAN of inner box | 4038917154611 |
| EAN case | 4038917027878 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

