



Product description

| | |
|---------------------|--|
| Order number | 94.6140.102 |
| Product description | x-well cell culture chamber, 1 well, on PCA slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister |

Product characteristics

| | |
|-----------------|-------------------|
| Closure type | slip-on lid |
| Design | on PCA slide |
| Cell type | Adherent cells |
| Number of wells | 1 |
| Growth area | 9 cm ² |
| Volume of work | 4 ml |
| Removable frame | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|-------|
| Width of product | 24 mm |
| Length of product | 76 mm |

Material & colours

| | |
|------------------|------------------------|
| Product material | Polycycloalkanes (PCA) |
| Frame material | Polystyrene (PS) |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Gas sterilisation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 6 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 6 |
| Piece (s) / outer case | 96 |
| Piece(s) / pallet | 10368 |
| Depth of box | 127 mm |
| Width of box | 19 mm |
| Height of box | 258 mm |
| Depth of case | 275 mm |
| Width of case | 270 mm |
| Height of case | 185 mm |
| Case volume | 0.0137 cbm |
| Weight of product | 0.0133 kg |
| Weight of case | 2.15 kg |
| EAN of inner box | 4038917303231 |
| EAN case | 4038917303224 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

