




Product Specification
x-well cell culture chamber, 1 well, on lumox®
slide, removable frame

Page 1



Product description

Order number	94.6150.101
Product description	x-well cell culture chamber, 1 well, on lumox® slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister

Product characteristics

Closure type	slip-on lid
Design	on lumox® slide
Cell type	Adherent cells
Number of wells	1
Growth area	9 cm ²
Volume of work	4 ml
Removable frame	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
x-well cell culture chamber, 1 well, on lumox®
slide, removable frame

Page 2

Size

Width of product	24 mm
Length of product	76 mm

Material & colours

Product material	lumox®
Frame material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	6
Type of smallest subpackaging	blister
Piece(s) / inner box	6
Piece (s) / outer case	96
Piece(s) / pallet	10368
Depth of box	127 mm
Width of box	19 mm
Height of box	258 mm
Depth of case	275 mm
Width of case	270 mm
Height of case	185 mm
Case volume	0.0137 cbm
Weight of product	0.0461 kg
Weight of case	5.3 kg
EAN of inner box	4038917255561
EAN case	4038917255554

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT