



Product description

Order number	83.1826.102
Product description	Syringe filter, Filtropur S plus, pre-filter: GF, membrane: CA, filtration surface: 6.2 cm ² , pore size: 0.2 µm, for sterile filtration, connection: Luer lock, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Example of use	for sterile filtration
Filter	yes
Surface area of filter	6.2 cm ²
Pore size	0,2
Connection	Luer lock

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Syringe filter, Filtropur S plus, CA, pore size: 0.2 µm, for sterile filtration

Page 2

Size

Diameter 33 mm

Material & colours

Product material Methacrylate-butadiene-styrene (MBS)

Colour of product light blue

Membrane material Cellulose acetate (CA)

Purity & certification

Product category not regulated

CE certified kein CE-Kennzeichen

Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic

Sterilisation Gamma irradiation

Pyrogen-free yes

Non-cytotoxic yes

Batched yes

Packaging

Minimum order qty. 50

Type of smallest subpackaging blister

Piece(s) / inner box 50

Piece (s) / outer case 500

Piece(s) / pallet 12000

Depth of box 200 mm

Width of box 133 mm

Height of box 107 mm

Depth of case 525 mm

Width of case 365 mm

Height of case 280 mm

Case volume 0.0537 cbm

Weight of product 0.0064 kg

Weight of case 4.2 kg

EAN of inner box 4038917133920

EAN case 4038917025263

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-03-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT