



Product description

Order number	83.1826.102
Product description	Syringe filter, Filtropur S plus, pre-filter: GF, membrane: CA, filtration surface: 6.2 cm ² , pore size: 0.2 µm, for sterile filtration, connection: Luer lock, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Example of use	for sterile filtration
Filter	yes
Surface area of filter	6.2 cm ²
Pore size	0,2
Connection	Luer lock

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Syringe filter, Filtropur S plus, CA, pore size: 0.2 µm, for sterile filtration

Page 2

Size

Diameter 33 mm

Material & colours

Product material Methacrylate-butadiene-styrene (MBS)
Colour of product light blue
Membrane material Cellulose acetate (CA)

Purity & certification

Product category no medical device | IVD
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation Gamma irradiation
Pyrogen-free yes
Non-cytotoxic yes
Batched yes

Packaging

Minimum order qty. 50
Type of smallest subpackaging blister
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 12000
Depth of case 375 mm
Width of case 285 mm
Height of case 530 mm
Case volume 0.0566 cbm
Weight of product 0.0064 kg
Weight of case 4.2 kg
EAN of inner box 4038917133920
EAN case 4038917025263

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

