



**Product Specification**  
**Filtropur BT 25, Bottle top filter, 250 ml, PES, 0.45  
µm**

**Page 1**



**Product description**

Order number	83.3940.100
Product description	Filtropur BT 25, Bottle top filter, for screw cap 45 mm, 250 ml, membrane: PES, Ø membraner: 90 mm, pore size: 0.45 µm, for clear filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

**Product characteristics**

Graduation	yes
Example of use	for clear filtration
Suitable for	screw cap 45 mm
Filter	yes
Surface area of filter	63.6 cm <sup>2</sup>
Pore size	0,45
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Sample volume	250 ml
Filtration volume	250 ml
Height of product	75.4 mm

### Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Membrane material	Polyethersulfone (PES)

### Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### Packaging

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	24
Piece(s) / pallet	576
Depth of case	521 mm
Width of case	203 mm
Height of case	381 mm
Case volume	0.0403 cbm
Weight of product	0.0752 kg
Weight of case	3.22 kg
EAN of inner box	4038917412179
EAN case	4038917412162

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

