



### **Product description**

Order number	83.3940.501
Product description	Filtropur V25, Vacuum filtration unit, 250 ml, screw cap, membrane: PES, Ø membraner: 50 mm, pore size: 0.22 µm, for sterile filtration, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

### **Product characteristics**

Graduation	yes
Closure type	screw cap
Example of use	for sterile filtration
Filter	yes
Surface area of filter	15.16 cm <sup>2</sup>
Pore size	0.22 µm
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	250 ml
Filtration volume	250 ml
Height of product	200 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Membrane material	Polyethersulfone (PES)

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	12
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	12
Piece(s) / pallet	420
Depth of case	433 mm
Width of case	358 mm
Height of case	236 mm
Case volume	0.0366 cbm
Weight of product	0.107 kg
Weight of case	2.72 kg
EAN of inner box	4038917969215
EAN case	4038917969222

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

