



### **Product description**

Order number	83.3940.505
Product description	Collection vessel, for Filtropur BT 25, 250 ml, screw cap, cap assembled, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

### **Product characteristics**

Graduation	yes
Closure type	screw cap
Suitable for	Filtropur BT 25
Volume of work	250 ml
Vessel type	flask
Filter	no

### **Size**

Sample volume	250 ml
Filtration volume	250 ml
Height of product	10.79 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	24
Piece(s) / pallet	840
Depth of case	433 mm
Width of case	358 mm
Height of case	236 mm
Case volume	0.0366 cbm
Weight of product	0.0552 kg
Weight of case	2.77 kg
EAN of inner box	4038917969369
EAN case	4038917969376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

