



## Product description

Order number	83.3940.511
Product description	Filtropur BT 25, Bottle top filter, for screw cap 45 mm, 250 ml, membrane: PES, Ø membraner: 50 mm, pore size: 0.22 µm, for sterile filtration, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

## Product characteristics

Graduation	yes
Example of use	for sterile filtration
Suitable for	screw cap 45 mm
Filter	yes
Surface area of filter	15.16 cm <sup>2</sup>
Pore size	0,22
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Filtropur BT 25, Bottle top filter, 250 ml, PES, 0.22  
µm**

**Page 2**

**Size**

Sample volume	250 ml
Filtration volume	250 ml
Height of product	11.35 cm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Membrane material	Polyethersulfone (PES)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	24
Piece(s) / pallet	840
Depth of case	433 mm
Width of case	358 mm
Height of case	236 mm
Case volume	0.0366 cbm
Weight of product	0.0565 kg
Weight of case	2.8 kg
EAN of inner box	4038917969406
EAN case	4038917969512

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

