



### Product description

Order number	83.3940
Product description	Filtropur V25, Vacuum filtration unit, 250 ml, screw cap, membrane: PES, Ø membraner: 90 mm, pore size: 0.45 µm, for clear filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

### Product characteristics

Graduation	yes
Closure type	screw cap
Example of use	for clear filtration
Filter	yes
Surface area of filter	63.6 cm <sup>2</sup>
Pore size	0,45
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Sample volume	250 ml
Filtration volume	250 ml
Height of product	143.9 mm

### Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polypropylene (PP)
Membrane material	Polyethersulfone (PES)

### Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### Packaging

Minimum order qty.	12
Type of inner packaging	bag
Piece(s) / inner box	1
Piece(s) / case	12
Piece(s) / pallet	288
Case volume	0.0405 cbm
Weight of product	0.0767 kg
Weight of case	2.72 kg
EAN of inner box	4038917411967
EAN case	4038917411950

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

