



Product description

Order number	83.3942.005
Product description	Collection vessel, for Filtropur BT 100, 1,000 ml, screw cap, cap assembled, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Suitable for	Filtropur BT 100
Volume of work	1.000 ml
Vessel type	flask
Filter	no

Size

Sample volume	1.000.00 ml
Filtration volume	1000 ml
Height of product	15.04 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polypropylene (PP)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	24
Type of inner packaging	bag
Piece(s) / inner box	1
Piece(s) / case	24
Piece(s) / pallet	288
Case volume	0.0863 cbm
Weight of product	0.1151 kg
Weight of case	5.31 kg
EAN of inner box	4038917412445
EAN case	4038917412438

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

