



### **Product description**

Order number	83.3942.005
Product description	Collection vessel, for Filtropur BT 100, 1,000 ml, screw cap, cap assembled, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

### **Product characteristics**

Graduation	yes
Closure type	screw cap
Suitable for	Filtropur BT 100
Volume of work	1.000 ml
Vessel type	flask
Filter	no

### **Size**

Sample volume	1.000.00 ml
Filtration volume	1000 ml
Height of product	15.04 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polypropylene (PP)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	24
Piece(s) / pallet	288
Depth of case	572 mm
Width of case	343 mm
Height of case	439 mm
Case volume	0.0861 cbm
Weight of product	0.1151 kg
Weight of case	5.31 kg
EAN of inner box	4038917412445
EAN case	4038917412438

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

