



Product description

Order number	83.3942
Product description	Filtropur V100, Vacuum filtration unit, 1,000 ml, screw cap, membrane: PES, Ø membraner: 90 mm, pore size: 0.45 µm, for clear filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Example of use	for clear filtration
Filter	yes
Surface area of filter	63.6 cm ²
Pore size	0,45
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	1.000.00 ml
Filtration volume	1000 ml
Height of product	28.33 cm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polypropylene (PP)
Membrane material	Polyethersulfone (PES)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	12
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	12
Piece(s) / pallet	144
Depth of case	572 mm
Width of case	343 mm
Height of case	439 mm
Case volume	0.0861 cbm
Weight of product	0.2593 kg
Weight of case	4.72 kg
EAN of inner box	4038917412117
EAN case	4038917412100

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

