



### **Product description**

Order number	83.3945.040
Product description	Cell strainer, pore size: 40 µm, blue, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

### **Product characteristics**

Example of use	for producing single-cell suspensions
Suitable for	50 ml centrifuge tubes
Filter	yes
Surface area of filter	4.2 cm <sup>2</sup>
Pore size	40

### **Size**

Diameter	34.09 mm
Height of product	21 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Colour of product	blue
Membrane material	nylon
Frame material	Polypropylene (PP)
Colour of frame	blue
Filter material	nylon

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	600
Piece(s) / pallet	9600
Case volume	0.0798 cbm
Weight of product	0.0068 kg
Weight of case	6.78 kg
EAN of inner box	4038917414982
EAN case	4038917414975

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

