



### **Product description**

Order number	83.3945.070
Product description	Cell strainer, pore size: 70 µm, white, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

### **Product characteristics**

Example of use	for producing single-cell suspensions
Suitable for	50 ml centrifuge tubes
Filter	yes
Surface area of filter	4.2 cm <sup>2</sup>
Pore size	70

### **Size**

Diameter	34.09 mm
Height of product	21 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Membrane material	nylon
Frame material	Polypropylene (PP)
Colour of frame	white
Filter material	nylon

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	50
Type of inner packaging	blister
Piece(s) / inner box	50
Piece(s) / case	600
Piece(s) / pallet	9600
Case volume	0.0798 cbm
Weight of product	0.0068 kg
Weight of case	6.78 kg
EAN of inner box	4038917415019
EAN case	4038917415002

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

