



Product description

Order number	83.3945.100
Product description	Cell strainer, pore size: 100 µm, yellow, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Example of use	for producing single-cell suspensions
Suitable for	50 ml centrifuge tubes
Filter	yes
Surface area of filter	4.2 cm ²
Pore size	100

Size

Diameter	34.09 mm
Height of product	21 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Membrane material	nylon
Frame material	Polypropylene (PP)
Colour of frame	yellow
Filter material	nylon

Purity & certification

Product category	not regulated
CE certified	kein CE-Kennzeichen
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	600
Piece(s) / pallet	9600
Depth of box	245 mm
Width of box	180 mm
Height of box	130 mm
Depth of case	520 mm
Width of case	385 mm
Height of case	410 mm
Case volume	0.0821 cbm
Weight of product	0.0068 kg
Weight of case	6.78 kg
EAN of inner box	4038917415040
EAN case	4038917415033

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

