




Product Specification
lumox® dish 50, Tissue culture dish, with foil
base, Ø: 50 mm, suspension cells

Page 1



Product description

Order number	94.6077.305
Product description	lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, for suspension cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Design	with foil base
Cell type	suspension cells
Base shape	flat base
Growth area	20.35 cm ²
Volume of work	5 ml - 10 ml
Vessel type	Dish
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**lumox® dish 50, Tissue culture dish, with foil
base, Ø: 50 mm, suspension cells**

Page 2

Size

Diameter	50 mm
Height of product	12 mm

Material & colours

Product material	lumox®
Frame material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	200
Piece(s) / pallet	7200
Depth of case	362 mm
Width of case	162 mm
Height of case	125 mm
Case volume	0.0073 cbm
Weight of product	0.01 kg
Weight of case	3.52 kg
EAN of inner box	4038917256070
EAN case	4038917256063

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

