



### Product description

|                     |   |
|---------------------|---|
| Order number        | 94.6000.034   |
| Product description | lumox® multiwell, Cell culture plate, with foil base, 384 well, frame: black, frame: PS, surface: adherent, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister, 4 piece(s)/inner box |

### Product characteristics

|                 |                      |
|-----------------|----------------------|
| Type of cap     | slip-on lid          |
| Design          | with foil base       |
| Cell type       | adherent             |
| Base shape      | flat base            |
| Number of wells | 384                  |
| Growth area     | 0,11 cm <sup>2</sup> |
| Volume of work  | 10 µl - 130 µl       |
| Tube            | Plate                |
| Coding colour   | red                  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

|                   |           |
|-------------------|-----------|
| Width of product  | 85,48 mm  |
| Height of product | 17 mm     |
| Length of product | 127,76 mm |

### Material & colours

|                  |                  |
|------------------|------------------|
| Product material | lumox®           |
| Frame material   | Polystyrene (PS) |
| Colour of frame  | black            |

### Purity & certification

|                  |  |
|------------------|--|
| Product category | no medical device   IVD                              |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Electron irradiation                                 |
| Pyrogen-free     | yes  |
| Batched          | yes  |

### Packaging

|                         |               |
|-------------------------|---------------|
| Minimum order qty.      | 4             |
| Type of inner packaging | blister       |
| Piece(s) / inner box    | 4             |
| Piece(s) / case         | 16            |
| Piece(s) / pallet       | 1024          |
| Depth of case           | 268 mm        |
| Width of case           | 191 mm        |
| Height of case          | 100 mm        |
| Case volume             | 0,0051 cbm    |
| Weight of product       | 0,0609 kg     |
| Weight of case          | 1,86 kg       |
| EAN of inner box        | 4038917254632 |
| EAN case                | 4038917254625 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

