



### **Product description**

Order number	94.6000.034
Product description	lumox® multiwell, Cell culture plate, with foil base, 384 well, frame: black, frame: PS, surface: adherent, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister, 4 piece(s)/inner box

### **Product characteristics**

Type of cap	slip-on lid
Design	with foil base
Cell type	adherent
Base shape	flat base
Number of wells	384
Growth area	0,11 cm <sup>2</sup>
Volume of work	10 µl - 130 µl
Tube	Plate
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Width of product	85,48 mm
Height of product	17 mm
Length of product	127,76 mm

### Material & colours

Product material	lumox®
Frame material	Polystyrene (PS)
Colour of frame	black

### Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

### Packaging

Minimum order qty.	4
Type of inner packaging	blister
Piece(s) / inner box	4
Piece(s) / case	16
Piece(s) / pallet	1024
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0,0051 cbm
Weight of product	0,0609 kg
Weight of case	1,86 kg
EAN of inner box	4038917254632
EAN case	4038917254625

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

