



Product description

Order number	94.6077.333
Product description	lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, for suspension cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Design	with foil base
Cell type	suspension cells
Base shape	flat base
Growth area	6.4 cm ²
Volume of work	2,5 ml
Vessel type	Dish
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, suspension cells

Page 2

Size

Diameter 35 mm
Height of product 6 mm

Material & colours

Product material lumox®
Frame material Polystyrene (PS)

Purity & certification

Product category no medical device | IVD
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation Electron irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 50
Type of inner packaging bag
Piece(s) / inner box 50
Piece(s) / case 250
Piece(s) / pallet 15000
Depth of case 348 mm
Width of case 253 mm
Height of case 265 mm
Case volume 0.0233 cbm
Weight of product 0.0037 kg
Weight of case 1.91 kg
EAN of inner box 4038917255141
EAN case 4038917255134

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-10-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT