



Product description

Order number	94.6077.410
Product description	lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, surface: adherent, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

Product characteristics

Type of cap	slip-on lid
Design	with foil base
Cell type	adherent
Base shape	flat base
Growth area	19,6 cm ²
Volume of work	5 ml - 10 ml
Tube	Dish
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, adherent

Page 2

Size

Diameter	50 mm
Height of product	12 mm

Material & colours

Product material	lumox®
Frame material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	50
Type of inner packaging	bag
Piece(s) / inner box	50
Piece(s) / case	200
Piece(s) / pallet	7200
Depth of case	362 mm
Width of case	162 mm
Height of case	125 mm
Case volume	0,0073 cbm
Weight of product	0,009 kg
Weight of case	4,59 kg
EAN of inner box	4038917255233
EAN case	4038917255226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-22

This document was prepared by EDP support and is valid without signature.



SARSTEDT