



Product description

Order number	94.6077.410
Product description	lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, for adherent cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Design	with foil base
Cell type	adherent cells
Base shape	flat base
Growth area	20.35 cm ²
Volume of work	5 ml - 10 ml
Vessel type	Dish
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
lumox[®] dish 50, Tissue culture dish, with foil base, Ø: 50 mm, adherent cells

Page 2

Size

Diameter 50 mm
Height of product 12 mm

Material & colours

Product material lumox[®]
Frame material Polystyrene (PS)

Purity & certification

Product category no medical device | IVD
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation Electron irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 50
Type of inner packaging bag
Piece(s) / inner box 50
Piece(s) / case 200
Piece(s) / pallet 7200
Depth of case 362 mm
Width of case 162 mm
Height of case 125 mm
Case volume 0.0073 cbm
Weight of product 0.009 kg
Weight of case 4.59 kg
EAN of inner box 4038917255233
EAN case 4038917255226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT