



### Product description

Order number	94.6110.024
Product description	lumox® multiwell, Cell culture plate, with foil base, 24 well, frame: black, frame: PS, surface: adherent, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister, 1 piece(s)/inner box

### Product characteristics

Type of cap	slip-on lid
Design	with foil base
Cell type	adherent
Base shape	flat base
Number of wells	24
Growth area	1,9 cm <sup>2</sup>
Volume of work	500 µl - 1500 µl
Tube	Plate
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**lumox® multiwell, Cell culture plate, with foil**  
**base, 24 well**

**Page 2**

**Size**

Width of product	85,48 mm
Height of product	23,4 mm
Length of product	127,76 mm

**Material & colours**

Product material	lumox®
Frame material	Polystyrene (PS)
Colour of frame	black

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	20
Type of inner packaging	blister
Piece(s) / inner box	1
Piece(s) / case	20
Piece(s) / pallet	1920
Depth of case	333 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0,0151 cbm
Weight of product	0,065 kg
Weight of case	1,67 kg
EAN of inner box	4038917255448
EAN case	4038917255431

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**