



Product description

| | |
|---------------------|--|
| Order number | 94.6077.138 |
| Product description | miniPERM® , Filling tube 5", for miniPERM® bioreactor, not autoclavable, sterile, 1 piece(s)/blister |

Product characteristics

| | |
|--------------|----------------------|
| Suitable for | miniPERM® bioreactor |
|--------------|----------------------|

Size

| | |
|-------------------|--------|
| Diameter | 5.9 mm |
| Width of product | 5.9 mm |
| Height of product | 127 mm |
| Length of product | 5.9 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Low Density Polyethylene (LD-PE) |
| Colour of product | transparent |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | sterile |
| Sterilisation | Gas sterilisation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 50 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 300 |
| Piece(s) / pallet | 14400 |
| Depth of case | 353 mm |
| Width of case | 268 mm |
| Height of case | 260 mm |
| Case volume | 0.0246 cbm |
| Weight of product | 0.001 kg |
| Weight of case | 1.55 kg |
| EAN of inner box | 4038917376488 |
| EAN case | 4038917255974 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

