



Product description

Order number 94.6001.054
Product description miniPERM® , Rack, for miniPERM® bioreactor, not autoclavable, 4 piece(s)/case

Product characteristics

Suitable for miniPERM® bioreactor

Material & colours

Product material Acrylnitrile butadiene styrene (ABS)
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	4
Type of smallest subpackaging	case
Piece(s) / inner box	4
Piece (s) / outer case	24
Piece(s) / pallet	672
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.085 kg
Weight of case	3.7 kg
EAN of inner box	4038917254786
EAN case	4038917254779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

