



### **Product description**

Order number 94.6001.061  
Product description Universal turning device, for 4 miniPERM® bioreactors, not autoclavable, 115/230 V, 1 piece(s)/case

### **Product characteristics**

Suitable for 4 miniPERM® bioreactors

### **Size**

Width of product 26 cm  
Length of product 46 cm

### **Purity & certification**

Product category no medical device | IVD  
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Packaging**

|                         |               |
|-------------------------|---------------|
| Minimum order qty.      | 1             |
| Type of inner packaging | case          |
| Piece(s) / inner box    | 1             |
| Piece(s) / case         | 1             |
| Piece(s) / pallet       | 30            |
| Depth of case           | 398 mm        |
| Width of case           | 358 mm        |
| Height of case          | 325 mm        |
| Case volume             | 0.0463 cbm    |
| Weight of product       | 5.861 kg      |
| Weight of case          | 6.3 kg        |
| EAN of inner box        | 4038917254878 |
| EAN case                | 4038917254861 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

