



### Product description

Order number 94.6001.094  
Product description miniPERM® , Start-up support kit, accessories for inoculation, sample collection and harvesting, for miniPERM® bioreactor, not autoclavable, set components: 8 x disposable syringe 50 ml, 20 x disposable syringe 2 ml, 8 x filling tube, 20 x needle, 6 x septum cap, 1 x miniPERM® stand, 1 piece(s)/bag

### Product characteristics

Application accessories for inoculation, sample collection and harvesting  
Suitable for miniPERM® bioreactor

### Purity & certification

Product category no medical device | IVD  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Packaging**

Minimum order qty.	1
Type of inner packaging	bag
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	96
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.481 kg
Weight of case	0.93 kg
EAN of inner box	4038917254939
EAN case	4038917254922

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

