



### Product description

Order number	94.6077.017
Product description	miniPERM® , HDC 50 production module, Tissue culture compartment, for protein, virus and biomass production, for suspension cells, volume of product module: 50 ml, not autoclavable, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

### Product characteristics

Closure type	screw cap
Design	for suspension cells
Application	for protein, virus and biomass production
Cell type	suspension cells

### Size

Production module volume	50 ml
Diameter	102 mm
Height of product	37 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polycarbonate (PC)
Closure material	Polypropylene (PP)
Dialysis membrane material	Cellulose acetate (CA)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	4
Type of inner packaging	bag
Piece(s) / inner box	4
Piece(s) / case	24
Piece(s) / pallet	480
Depth of case	498 mm
Width of case	364 mm
Height of case	310 mm
Case volume	0.0562 cbm
Weight of product	0.1405 kg
Weight of case	5.11 kg
EAN of inner box	4038917385626
EAN case	4038917385602

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

