



Product description

Order number	94.6077.135
Product description	Luer cannula, 25G x 1", coding colour: orange, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 100 piece(s)/case

Product characteristics

Coding colour	orange
Safety product	no
Connection	Luer (f)

Size

Diameter of needle (external)	25 G / 0,5 mm
Length of needle	25 mm

Purity & certification

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Batched yes

Packaging

Minimum order qty.	100
Type of inner packaging	blister
Piece(s) / inner box	100
Piece(s) / case	1200
Piece(s) / pallet	129600
Depth of case	383 mm
Width of case	193 mm
Height of case	190 mm
Case volume	0.014 cbm
Weight of product	0.0011 kg
Weight of case	1.8 kg
EAN of inner box	4038917255899
EAN case	4038917255882

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

