



### **Product description**

Order number	94.6150.801
Product description	x-well cell culture chamber, 8-well, on lumox® slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister

### **Product characteristics**

Closure type	slip-on lid
Design	on lumox® slide
Cell type	adherent cells
Number of wells	8
Growth area	0.8 cm <sup>2</sup>
Volume of work	0,5 ml
Removable frame	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Width of product	24 mm
Length of product	76 mm

### Material & colours

Product material	lumox®
Frame material	Polystyrene (PS)

### Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

### Packaging

Minimum order qty.	6
Type of smallest subpackaging	blister
Piece(s) / inner box	6
Piece (s) / outer case	96
Piece(s) / pallet	10368
Depth of box	127 mm
Width of box	19 mm
Height of box	258 mm
Depth of case	275 mm
Width of case	270 mm
Height of case	185 mm
Case volume	0.0137 cbm
Weight of product	0.0479 kg
Weight of case	5.48 kg
EAN of inner box	4038917256223
EAN case	4038917256216

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

