



Product description

Order number 90.183.100
Product description Tube holder, SC 2700, for centrifuge SC 2700

Product characteristics

Suitable for centrifuge SC 2700

Size

Diameter 18.4 mm
Width of product 23 mm
Height of product 83.5 mm
Length of product 23 mm

Material & colours

Colour of product black

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2
Type of smallest subpackaging bag
Piece(s) / inner box 2
Piece (s) / outer case 50
Piece(s) / pallet 16200
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0.0051 cbm
Weight of product 0.021 kg
Weight of case 1.2 kg
EAN of inner box 4038917909143
EAN case 4038917909150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

