



Product description

Order number 92.184.555
Product description Cap, LC 6, natural, for Tube holder

Product characteristics

Suitable for Tube holder

Material & colours

Colour of product natural

Purity & certification

Product in accordance with MDD or IVDD no medical device | IVD
Lot or serial number no

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 60
Piece(s) / pallet 19440
Depth of case 268 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Cap, LC 6, natural, for Tube holder

Page 2

Width of case	191 mm
Height of case	100 mm
Case volume	0,0051 cbm
Weight of product	0,0039 kg
Weight of case	0,43 kg
EAN of inner box	4038917102483
EAN case	4038917073257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-06-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT