



### **Product description**

Order number 92.184.755  
Product description Cap, LC 24, natural, for Tube holder

### **Product characteristics**

Suitable for Tube holder

### **Material & colours**

Colour of product natural

### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 1  
Piece(s) / inner box 1  
Piece(s) / case 1  
Piece(s) / pallet 324  
Depth of case 268 mm  
Width of case 191 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Cap, LC 24, natural, for Tube holder**

**Page 2**

Height of case	100 mm
Case volume	0,0051 cbm
Weight of product	0,076 kg
Weight of case	0,21 kg
EAN of inner box	4038917102568
EAN case	4038917073318

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**