



Product Specification

Centrifuge, SC 2700, supply voltage: 115 - 230 V,
2 preset programmes and 1 variable programme

Page 1



Product description

Order number	90.183.000
Product description	Centrifuge, SC 2700, supply voltage: 115 - 230 V, centrifugation up to max. 2,700 x g, incl. swing-out rotor 6-fold and tube holder 1-fold for tubes up to 100 mm length, for separation of serum/plasma from the blood cells, 2 preset programmes and 1 variable programme

Product characteristics

Design	2 preset programmes and 1 variable programme
Suitable for	separation of serum/plasma from the blood cells
Centrifugation max (RCF)	2700 x g
Number of positions	6

Size

Width of product	350 mm
Height of product	260 mm
Length of product	430 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Centrifuge, SC 2700, supply voltage: 115 - 230 V,
2 preset programmes and 1 variable programme**

Page 2

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	20
Depth of box	540 mm
Width of box	400 mm
Height of box	385 mm
Depth of case	540 mm
Width of case	400 mm
Height of case	385 mm
Case volume	0.0772 cbm
Weight of product	16.655 kg
Weight of case	18.5 kg
EAN of inner box	4038917574877
EAN case	4038917574860

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT