



Product description

Order number	90.183.000
Product description	Centrifuge, SC 2700, supply voltage: 115/230 V, centrifugation up to max. 2700 x g, incl. swing-out rotor 6-fold and tube holder 1-fold for tubes up to x 100 mm length, for separation of serum/plasma from the blood cells

Product characteristics

Design	2 voreingestellte Programme und 1 variables Programm
Suitable for	separation of serum/plasma from the blood cells
Centrifugation max (RCF)	2700 x g
Number of positions	6

Size

Width of product	350 mm
Height of product	260 mm
Length of product	430 mm

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE supplier
Batched	no

Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	20
Case volume	0,0772 cbm
Weight of product	16,655 kg
Weight of case	18,5 kg
EAN of inner box	4038917574877
EAN case	4038917574860

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

