



**Product Specification**  
**length: 1.5 m, transparent**

**Page 1**

**Product description**

Order number 90.2620.301  
Product description length: 1.5 m, transparent, 1 piece(s)/bag

**Product characteristics**

Suitable for IVARO

**Size**

Diameter 0.8 mm  
Length of product 1.5 m

**Material & colours**

Product material Bioprene  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging bag  
Piece(s) / inner box 1  
Piece (s) / outer case 10  
Case volume 0.0062 cbm  
Weight of product 0.0005 kg  
EAN of inner box 4038917584234  
EAN case 4038917584227

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**