



Product description

Order number 90.2620.301
Product description length: 1.5 m, transparent, 1 piece(s)/bag

Product characteristics

Suitable for IVARO

Size

Diameter 0.8 mm
Length of product 1.5 m

Material & colours

Product material Bioprene
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging bag
Piece(s) / inner box 1
Piece (s) / outer case 10
Case volume 0.0062 cbm
Weight of product 0.0005 kg
EAN of inner box 4038917584234
EAN case 4038917584227

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

