



Product description

Order number 92.170.153
Product description Wash bottle, 100 ml

Size

Sample volume 100 ml
Max. volume 100 ml

Material & colours

Product material borosilicate glass
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Case volume 0,0007 cbm
Weight of product 0,002 kg
Weight of case 0,05 kg
EAN of inner box 4038917100953
EAN case 4038917051293

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

