



Product Specification

Sample changer B/G type GS 301, for gas sampler GS 301 (item no 90.170.350) & tubes 7 x 125 mm

Page 1



Product description

Order number	92.170.352
Product description	Sample changer B/G type GS 301, for gas sampler GS 301 (item no 90.170.350) & tubes 7 x 125 mm, 1 piece(s)/case

Product characteristics

Suitable for	gas sampler GS 301 (item no 90.170.350) & tubes 7 x 125 mm
--------------	--

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0.0011 cbm
Weight of product	0.6 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-12-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Sample changer B/G type GS 301, for gas sampler
GS 301 (item no 90.170.350) & tubes 7 x 125 mm**

Page 2

Weight of case	0.7 kg
EAN of inner box	4038917101196
EAN case	4038917051514

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-12-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT