



Product description

Order number	92.170.352
Product description	Sample changer B/G type GS 301, for tubes 7 x 125 mm, 1 piece(s)/inner box

Product characteristics

Suitable for	tubes 7 x 125 mm
--------------	------------------

Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

Packaging

Minimum order qty.	1
Type of inner packaging	inner box
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0011 cbm
Weight of product	0,6 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Sample changer B/G type GS 301, for tubes 7 x 125 mm

Weight of case	0,7 kg
EAN of inner box	4038917101196
EAN case	4038917051514

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

