



### Product description

Order number	92.170.353
Product description	Sample changer N type GS 301, for tubes 6 x 70 mm, 1 piece(s)/case

### Product characteristics

Suitable for	tubes 6 x 70 mm
--------------	-----------------

### Purity & certification

Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0.0009 cbm
Weight of product	0.5 kg
Weight of case	0.65 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Sample changer N type GS 301, for tubes 6 x 70 mm**

**Page 2**

EAN of inner box                      4038917101202  
EAN case                                 4038917051521

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-08-15

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**