



Product description

Order number	92.170.354
Product description	Sample changer ATD type GS 301, for tubes 6 x 89 mm, 1 piece(s)/inner box

Product characteristics

Suitable for	tubes 6 x 89 mm
--------------	-----------------

Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

Packaging

Minimum order qty.	1
Type of inner packaging	inner box
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0011 cbm
Weight of product	0,55 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Sample changer ATD type GS 301, for tubes 6 x 89 mm

Page 2

Weight of case	0,7 kg
EAN of inner box	4038917101219
EAN case	4038917051538

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-05-10

This document was prepared by EDP support and is valid without signature.



SARSTEDT